

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

**Product Name:**

Insulin-like Growth Factor-I from mouse - recombinant, expressed in Escherichia coli, lyophilized powder, cell culture tested

**Product Number:** 18779**MDL:** MFCD00081716**Storage Temperature:** -20 °C

TEST	Specification
Appearance (Form) Lyophilized from a 0.2 µm-filtered solution in PBS	Lyophilized Powder
Bioactivity EC50 (Cell proliferation, MCF-7 cells) Report Result	
Purity > 97% by SDS-PAGE	Conforms
Endotoxin Assay < 1.0 EU/µg	Conforms

Specification: PRD.0.ZQ5.10000012820

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.