


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Methyl linolenate - ≥99% (GC)

**Product Number:** L2626  
**CAS Number:** 301-00-8  
**MDL:** MFCD00135851  
**Formula:** C<sub>19</sub>H<sub>32</sub>O<sub>2</sub>  
**Formula Weight:** 292.46 g/mol  
**Storage Temperature:** -20 °C



### TEST

### Specification

Appearance (Color)	Colorless to Very Faint Yellow
Appearance (Turbidity)	Clear
Appearance (Form)	Liquid
<sup>13</sup> C NMR Identity	Conforms to Structure
Proton NMR Spectrum	Conforms to Structure
Purity (GC)	≥ 99 %
Recommended Retest Period	-----
6 years	

Specification: PRD.2.ZQ5.10000035881

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.