

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Heptanone - 99%

Product Number: 537683
CAS Number: 110-43-0
MDL: MFCD00009513
Formula: C₇H₁₄O
Formula Weight: 114.19 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Refractive index at 25 ° C	1.406 - 1.409
Water (by Karl Fischer)	≤ 0.05 %
Specific Gravity at 20 Degrees Celsius	0.815 - 0.818
Acidity as Acetic Acid	≤ 0.01 %
Assay	≥ 99.0 %
Impurity % Other Ketones	≤ 0.50 %
Alcohol Content	≤ 0.20 %

Specification: PRD.0.ZQ5.10000030920

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.