

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Potassium hydroxide on alumina

Product Number: 720143

KOH

Formula Weight: 56.11 g/mol

TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Powder
ICP Major Analysis Confirms Potassium and Aluminum Components	Confirmed
Titration with ZnSO ₄ % Al Specification: Report Result	
Loss on Drying	≤ 5.0 %

Specification: PRD.0.ZQ5.10000074046

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.