

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Calcium bromide hydrate – Vetec™ reagent grade, 98%

Product Number:

V900070

CAS Number:

71626-99-8



Formula:

 $\text{Br}_2\text{Ca} \cdot x\text{H}_2\text{O}$

Formula Weight:

199.90 g/mol

TEST

Specification

Appearance (Colour)

White

Appearance (Form)

Powder or Chunks

ICP Major Analysis

Confirmed

Confirms Calcium Component

Titration by AgNO₃

97.5 - 102.5 %

As xH₂O

Specification: PRE.1.ZQ5.10000002585

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.