


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Graphene oxide - 2 mg/mL, dispersion in H₂O

Product Number: 763705



TEST	Specification
------	---------------

Appearance (Color)	Brown to Very Dark Brown and Black
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Carbon (Dry Basis)	42.0 - 52.0 %
Sulfur (Dry Basis)	≤ 4.0 %
Potassium (K)	≤ 60 ppm
Manganese (Mn)	≤ 60 ppm
pH	
Specification: Report Only	
Residue on Evaporation	1.9 - 2.4
mg/ml	

Specification: PRD.3.ZQ5.10000030262

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.