

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

2',7'-Dichlorofluorescein diacetate - ≥97%

Product Number:

D6883

CAS Number:

4091-99-0

Formula:


C<sub>24</sub>H<sub>16</sub>Cl<sub>2</sub>O<sub>7</sub>

Formula Weight:

487.29 g/mol

Storage Temperature:

-20 °C



### TEST

### Specification

Appearance (Form)	Powder
Appearance (Colour)	White to White w/ Yellow Cast
Solubility (Solvent)	Methanol
Solubility (Conc)	9.80 - 10.20 mg/ml
Solubility (Heating)	Yes or No
Solubility (Turbidity)	Clear
Solubility (Color)	Colorless to Faint Yellow
Elemental Anal. (%C anhydrous)	58.65 - 59.65 %
Purity (TLC)	≥ 97.00 %
NMR (Solvent)	DMSO-d <sub>6</sub>
Identity by NMR	Consistent
Purity (HPLC)	≥ 97.00 %
Storage Conditions	Under Argon
Recommended Retest Period	-----
3 Years	

Specification: PRE.3.ZQ5.10000001064

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.