

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Hyperbranched bis-MPA polyester-64-hydroxyl, generation 4 - ≥97%

Product Number: 686573
CAS Number: 326794-48-3
Formula: C315H512O189
Formula Weight: 7,323.32 AMU

TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Carbon	49.1 - 52.1 %
Molecular Weight by GPC Specification: Report Results	
Molecular Number by GPC Specification: Report Results	
Poly Dispersity MW/MN Specification: Report Results	
Note May Contain Trace Amounts of Methanol	Note

Specification: PRD.0.ZQ5.10000029112

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.