

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tin(II) sulfide - 96%

Product Number: 367168
CAS Number: 1314-95-0
MDL: MFCD00011245
Formula: SSn
Formula Weight: 150.78 g/mol

SnS

TEST	Specification
Appearance (Color) Grey	Conforms to Requirements
Appearance (Form) Crystalline Powder, Flakes, Chips or Chunks	Conforms to Requirements
X-Ray Diffraction	Conforms to Structure
Complexometric EDTA % Sn	75.2 - 82.3 %

Specification: PRD.0.ZQ5.10000040584

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.