

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Tin(IV) oxide - nanopowder,  $\leq 100$  nm avg. part. size**Product Number:****549657**

CAS Number:

18282-10-5

MDL:

MFCD00011244

Formula:

O<sub>2</sub>Sn

Formula Weight:

150.71 g/mol



TEST	Specification
------	---------------

Appearance (Color)

White to Grey

Appearance (Form)

Powder

X-Ray Diffraction

Conforms

Surface Area

Specification: Report Results

Average Particle Size

 $\leq 100$  nm

Product of Supplier

Confirmed

Engi-Mat Co.

Specification: PRD.4.ZQ5.10000021202

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.