

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

## **Product Specification**

Product Name:

lodine - ACS reagent, ≥99.8%, solid

 Product Number:
 207772

 CAS Number:
 7553-56-2

 MDL:
 MFCD00011355

Formula: 12

Formula Weight: 253.81 g/mol

12

## TEST Specification

Appearance (Color)
Appearance (Form)
Titration by Na2S2O3
Nonvolatile Matter

Chlorine and Bromine (as Cl) Meets ACS Requirements Grey to Very Dark Grey Beads or Flakes ≥ 99.8 % ≤ 0.01 %

< = 0.005%

Current ACS Specification

Specification: PRD.6.ZQ5.10000003931

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.