


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Product Number: 393304
CAS Number: 29823-18-5
MDL: MFCD00037103
Formula: C₉H₁₇BrO₂
Formula Weight: 237.13 g/mol



TEST	Specification
Appearance (Color)	Colorless to Yellow
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Titration by AgNO ₃ % Purity Based On Bromide Content	96.5 - 103.5 %
By Titration With AgNO ₃ After O ₂ Combustion	
Purity (GC)	≥ 96.5 %

Specification: PRD.0.ZQ5.10000082721

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.