

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
2,8,9-Triisopropyl-2,5,8,9-tetraaza-1-phosphabicyclo[3,3,3]undecane

Product Number: 556955
CAS Number: 175845-21-3
MDL: MFCD03701530
Formula: C₁₅H₃₃N₄P
Formula Weight: 300.42 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
White to Yellow (Moist Crystalline Solid) or Colorless to Yellow (Liquid)	
Appearance (Form)	Conforms to Requirements
Moist Crystalline Solid or Liquid	
Infrared spectrum	Conforms to Structure
Proton NMR Spectrum	Conforms to Structure
Phosphorus NMR Spectra	Conforms to Structure

Specification: PRD.1.ZQ5.10000054015

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.