


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2,4-Dichloro-5-fluoropyrimidine - 97%

Product Number: **653233**
CAS Number: 2927-71-1
MDL: MFCD00233551
Formula: C₄HCl₂FN₂
Formula Weight: 166.97 g/mol



| TEST | Specification |
|------|---------------|
|------|---------------|

| | |
|---|--|
| Appearance (Color) | White to Light Yellow and Light Beige to Beige |
| Appearance (Form) | Powder or Crystals or Solid or Chunks |
| Infrared spectrum | Conforms to Structure |
| Titration with AgNO ₃ After Oxygen Combustion | 96.5 - 103.5 % |
| Purity (GC) | ≥ 96.5 % |

Specification: PRD.3.ZQ5.10000018085

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

