

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Calcium oxalate - ≥99.9% trace metals basis

**Product Number:** 455997  
**CAS Number:** 563-72-4  
**MDL:** MFCD00012474  
**Formula:** C<sub>2</sub>CaO<sub>4</sub>  
**Formula Weight:** 128.10 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Chunks
X-Ray Diffraction	Conforms to Structure
Purity	Confirmed
> = 99.9% Purity Based On Trace Metals Analysis	
Trace Metal Analysis	≤ 1000 ppm

Specification: PRD.2.ZQ5.10000066526

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.