

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Anti-TRPC6 antibody produced in rabbit – affinity isolated antibody, lyophilized powder

Product Number:

T6442

MDL:

MFCD02686669

TEST	Specification
------	---------------

Formulation	Conforms
Lyophilized from phosphate buffered saline, pH 7.4, containing 1% BSA, and 0.05% sodium azide	
Concentration	Conforms
> or = 0.5 mg/mL after reconstitution	
Working Dilution	Conforms
by Direct Immunoblotting	
1:200 to 1:400	
using rat brain lysates, HRP goat anti-rabbit and chemiluminescent detection	
Immunohistochemistry	Conforms
suitable using mouse brain sections	

Specification: PRD.1.ZQ5.10000009485

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.