

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Trypsin from human pancreas – salt-free, lyophilized powder, vial of $\geq 1,000$ BAEE units

Product Number:

T6424

CAS Number:

9002-07-7

MDL:

MFCD01323069

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

 μg Protein/Vial ≥ 75

units/vial

 ≥ 1000

One BAEE unit will produce a delta A253 of 0.001 per minute in a reaction volume of 3.2 mL at pH 7.6 at 25 deg C. (1 cm light path)

BTEE units/vial

 ≤ 0.10

Chymotrypsin

Tested For Infectious Agents

Tested

Specification: PRD.1.ZQ5.10000013467

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.