


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
cis-Bis(benzonitrile)dichloroplatinum(II) - 98%

Product Number: 275816
CAS Number: 15617-19-3
MDL: MFCD00013125
Formula: C₁₄H₁₀Cl₂N₂Pt
Formula Weight: 472.23 g/mol



TEST

Specification

Appearance (Color)	Yellow
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Carbon	34.7 - 36.5 %
Nitrogen	5.7 - 6.1 %
ICP Major Analysis	Confirmed
Confirms Platinum Component	

Specification: PRD.0.ZQ5.10000038281

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.