

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

N,N'-Diphenyl-p-phenylenediamine - 98%

Product Number:**292265**

CAS Number:

74-31-7

MDL:


MFCD00003015

Formula:

C₁₈H₁₆N₂

Formula Weight:

260.33 g/mol

**TEST****Specification**

Appearance (Color)

Faint Grey to Very Dark Grey and Brown

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Carbon

81.0 - 85.1 %

Nitrogen

10.4 - 11.1 %

Specification: PRD.0.ZQ5.10000049029

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.