

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Phytamax™ Orchid Maintenance Medium - powder, plant cell culture tested

Product Number: **P6668**

Storage Temperature: 2 - 8 °C

TEST	Specification
Solubility	Conforms
Partially Soluble in water	
pH	4.75 - 5.75
27.31 g/L	
Appearance (Color)	Conforms
Under Review	
Appearance (Form)	Conforms
Fine to Coarse Powder	
Solubility (Turbidity)	Conforms
Opaque, Grainy	
Solubility (Color)	Conforms
Black	
Plant Cell Culture Test	Pass
Supports and/or facilitates plant growth and/or shoot proliferation in two or more plant tissue cultured lines with no morphological aberrations to plants.	

Specification: PRD.0.ZQ5.10000048577

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.