

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

N,O-Bis(trimethylsilyl)acetamide - synthesis grade, ≥95%

 Product Number:
 128910

 CAS Number:
 10416-59-8

 MDL:
 MFCD00008270

 Formula:
 C8H21NOSi2

 Formula Weight:
 203.43 g/mol

TEST Specification

Appearance (Color) Appearance (Form) Infrared spectrum Purity (GC)

Specification: PRD.1.ZQ5.10000021551

Colorless to Light Yellow Liquid Conforms to Structure > 95.0 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.