


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Triphenyltin hydride

Product Number: 244775
CAS Number: 892-20-6
MDL: MFCD00003004
Formula: C₁₈H₁₆Sn
Formula Weight: 351.03 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)
Appearance (Form)
Proton NMR Spectrum

Colorless to Yellow or White to Yellow
Solid or Liquid
Conforms to Structure

Specification: PRD.0.ZQ5.10000050905

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.