

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Hexane - ReagentPlus*, ≥99%

Product Number: 139386 CAS Number: 110-54-3

MDL: MFCD00009520

Formula: C6H14
Formula Weight: 86.18 g/mol

CH₃(CH₂)₄CH₃

TEST

Appearance (Color)

Appearance (Form) Infrared spectrum Refractive index at 20 ° C

Purity (GC)

Specification: PRD.0.ZQ5.10000028065

Specification

Colorless

Liquid
Conforms to Structure
1.373 - 1.377
≥ 99.0 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.