

Product Specification

Product Name:

Gentamicin sulfate salt - powder, BioReagent, suitable for cell culture

Product Number:

G1264

CAS Number:

1405-41-0



Storage Temperature:

2 - 8 °C

TEST

Specification

Appearance (Form)

Powder

Appearance (Color)

White to Off White

Water (by Karl Fischer)

≤ 15 %

Assay (Anhydrous Basis)

≥ 590 IU/mg

Cell Culture Test

Pass

Expiration Date Period

4 Years

Specification: PRE.1.ZQ5.10000020069

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

