

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Protein Kinase A inhibitor fragment 14-22, myristoylated trifluoroacetate salt - ≥90% (HPLC), lyophilized powder

Product Number:

P9115

MDL:


MFCD03788010

Formula:

C₅₃H₉₈N₂₀O₁₂ · C₂HF₃O₂

Formula Weight:

1,321.50 g/mol



| TEST | Specification |
|---------------------------|---------------|
| Appearance (Color) | White |
| Appearance (Form) | Powder |
| Solubility (Color) | Colorless |
| Solubility (Turbidity) | Clear |
| 1 mg/ml, H ₂ O | |
| Purity (HPLC) | ≥ 90 % |

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.