

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium persulfate - ACS reagent, ≥99.0%

 Product Number:
 216224

 CAS Number:
 7727-21-1

 MDL:
 MFCD00011386

Formula: K2S2O8 Formula Weight: 270.32 g/mol

TEST Specification

Appearance (Color) White

Appearance (Form) Powder or Crystals
X-Ray Diffraction Conforms to Structure

Titration with KMnO4 \geq 99.0 %

Insoluble Matter < 0.005 % c = 10%; Water

Chlorine Compounds < 0.001 %

 Iron (Fe)
 ≤ 5 ppm

 Heavy Metal (as Lead)
 ≤ 0.001 %

Manganese < 2 ppm

Meets ACS Requirements Current ACS Specification

Specification: PRD.1.ZQ5.10000042090

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.