

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Hyaluronic acid sodium salt from bovine vitreous humor

Product Number: **H7630**
CAS Number: 9067-32-7
MDL: MFCD00131348

Storage Temperature: -20 °C


TEST

Appearance (Color)
Appearance (Form)
Solubility (Color)
Solubility (Turbidity)
0.3 mg/ml, 0.3 M NaH₂PO₄ buffer, pH 5.35
Sodium (Na)
Potassium (K)
Suitability
Suitable as a hyaluronidase substrate

Specifi

White
Powder
Colorless
Clear

3 - 10 %
≤ 1 %
Suitable



Specification: PRD.0.ZQ5.10000039263

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.