

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Glycopeptidase A from almonds - buffered aqueous glycerol solution,  $\geq 0.05$  unit/mL

Product Number:

**G0535**

CAS Number:

83534-39-8

MDL:

MFCD00131222

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

units/ml	$\geq 0.05$
Enzymatic Activity	
One unit will hydrolyze 1.0 micromole of ovalbumin glycopeptide per minute at pH 5.0 at 37°C.	
Formulation	Conforms
Solution in 50% glycerol containing 50 mM citrate-phosphate buffer, pH 5.0, containing BSA.	
ug BSA/mL	$\geq 20$

Specification: PRD.1.ZQ5.10000012251

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.