


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Polypropylene – Isotactic, average M_w ~340,000, average M_n ~97,000

Product Number: 427861
CAS Number: 9003-07-0
MDL: MFCD00084447
Formula: C₃H₆



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Beads
Infrared spectrum	Conforms to Structure
Melt Index (g/10 min 230C/2.16Kg)	3.0 - 5.0

Specification: PRD.0.ZQ5.10000101198

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.