

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Product Number: 208086
CAS Number: 7646-85-7
MDL: MFCD00011295
Formula: Cl₂Zn
Formula Weight: 136.30 g/mol

ZnCl₂

TEST	Specification
------	---------------

Appearance (Color)	White to Off-White
Appearance (Form)	Powder or Crystal or Chunk(s)
X-Ray Diffraction	Conforms to Structure
Compleximetric EDTA % Zn	≥ 46.9 %

Specification: PRD.0.ZQ5.10000032223

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.