

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium phosphate dibasic trihydrate - ReagentPlus®, ≥99.0%

Product Number:

P5504



CAS Number:

16788-57-1

MDL:

MFCD00149926

Formula:

HK2O4P · 3H2O

Formula Weight:

228.22 g/mol

TEST**Specification**

Appearance (Color)

Colorless or White

Appearance (Form)

Powder or Crystals

Solubility (Color)

Colorless

Solubility (Turbidity)

Clear

100 mg/mL, H₂O

Titration with HCl

≥ 99.0 %

Recommended Retest Period

3 years

Specification: PRD.1.ZQ5.10000015785

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.