


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Ethylene glycol – anhydrous, 99.8%

Product Number: 324558
CAS Number: 107-21-1
MDL: MFCD00002885
Formula: C₂H₆O₂
Formula Weight: 62.07 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 99.75 %
Water (by Karl Fischer)	≤ 0.003 %
Note	-----

For 100 mL unit size, water content may exceed above specification; not to exceed 0.005%.

Specification: PRD.2.ZQ5.10000004111

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.