


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Diiodomethane - ReagentPlus®, 99%, contains copper as stabilizer

Product Number: 158429
CAS Number: 75-11-6
MDL: MFCD00001079
Formula: CH₂I₂
Formula Weight: 267.84 g/mol



| TEST | Specification |
|------|---------------|
|------|---------------|

| | |
|--------------------|-------------------------------------|
| Appearance (Color) | Colorless to Green-Yellow to Yellow |
| Appearance (Form) | Liquid |
| Infrared Spectrum | Conforms to Structure |
| Purity (GC) | ≥ 98.5 % |

Specification: PRD.2.ZQ5.10000018900

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.