


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1,2-Diaminopropane-N,N,N',N'-tetraacetic acid - ≥98%

Product Number: 158135
CAS Number: 4408-81-5
MDL: MFCD00004286
Formula: C₁₁H₁₈N₂O₈
Formula Weight: 306.27 g/mol



| TEST | Specification |
|--------------------|-----------------------|
| Appearance (Color) | White |
| Appearance (Form) | Powder |
| Infrared spectrum | Conforms to Structure |
| Carbon | 42.3 - 44.0 % |
| Nitrogen | 8.8 - 9.5 % |

Specification: PRD.1.ZQ5.10000042286

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.