

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-phospho-Insulin Receptor (IR) (pTyr972) antibody produced in rabbit - affinity isolated antibody, buffered

aqueous solution

Product Number: I1783

MDL: MFCD05665246

Storage Temperature: -70 ℃

TEST Specification

Formulation Conforms

 $Solution in \ Dulbecco's \ Phosphate \ Buffered \ Saline \ (without \ Mg2+ \ and \ Ca2+), \ pH \ 7.3, \ 50\% \ Glycerol \ with \ 1 \ mg/ml \ BSA \ and \ 0.05\%$

Sodium Azide,

Western Blot Conforms

assay dependent Sufficient for 10 blots

Specification: PRD.1.ZQ5.10000067798

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.