


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Boron trifluoride tetrahydrofuran complex

Product Number: 434280
CAS Number: 462-34-0
MDL: MFCD00040372
Formula: C₄H₈BF₃O
Formula Weight: 139.91 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Colorless to Yellow to Brown
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Purity (Titration by NaOH) %BF ₃	47.5 - 49.5 %

Specification: PRD.0.ZQ5.10000068396

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.