

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Dichlorobis(tri-*o*-tolylphosphine)palladium(II) - 97%

Product Number:

453048

CAS Number:

40691-33-6

MDL:


MFCD00274659

Formula:

C₄₂H₄₂Cl₂P₂Pd

Formula Weight:

786.06 g/mol



TEST

Specification

Appearance (Color)

Yellow to Orange

Appearance (Form)

Powder or Crystals

Infrared spectrum

Conforms to Structure

Carbon

61.9 - 66.4 %

Specification: PRD.0.ZQ5.10000021208

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.