

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

4-Methoxy-2-methylphenylboronic acid

Product Number:

639370

CAS Number:

208399-66-0

MDL:


MFCD02684315

Formula:

C₈H₁₁BO₃

Formula Weight:

165.98 g/mol

**TEST****Specification**

Appearance (Color)

White to Off White

Appearance (Form)

Powder

Proton NMR spectrum

Conforms to Structure

Purity (Titration by NaOH)

≥ 95.0 %

Specification: PRD.0.ZQ5.10000037917

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.