


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Amino-3-methylpyridine - 95%

Product Number: A75633
CAS Number: 1603-40-3
MDL: MFCD00006320
Formula: C₆H₈N₂
Formula Weight: 108.14 g/mol



TEST	Specification
Appearance (Color) Yellow to Dark Yellow	Conforms to Requirements
Appearance (Form) Moist Crystals or Solid or Liquid	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 94.5 %

Specification: PRD.1.ZQ5.10000039098

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.