

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Anti-TRPC3 antibody produced in rabbit – affinity isolated antibody, lyophilized powder

Product Number: **T5067**

MDL: MFCD02686005

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Formulation	Conforms
-------------	----------

Lyophilized from PBS, pH 7.4,  
containing 1% BSA and 0.05% sodium azide

Concentration	Conforms
---------------	----------

> /= 0.6 mg/ml

Working Dilution	Conforms
------------------	----------

by direct Immunoblotting  
1:200 to 1:1000

Specification: PRD.0.ZQ5.10000100140

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.