


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Phenylphenol - 85%

Product Number: 262250
CAS Number: 580-51-8
MDL: MFCD00002294
Formula: C₁₂H₁₀O
Formula Weight: 170.21 g/mol

**TEST****Specification**

Appearance (Color)	White to Beige
Appearance (Form)	Powder or Crystals or Chunk(s) or Powder with Chunk(s) or Crystals with Chunk(s)
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 83.0 %

Specification: PRD.1.ZQ5.10000018106

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.