

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Copper(I) iodide – 99.999% trace metals basis

Product Number: 215554
CAS Number: 7681-65-4
MDL: MFCD00010978
Formula: CuI
Formula Weight: 190.45 g/mol

CuI

TEST	Specification
Appearance (Color)	Brown-Grey
Appearance (Form)	Powder or Chunks
Complexometric EDTA % Copper	31.0 - 34.0 %
ICP Major Analysis Confirms Copper Component	Confirmed
Purity 99.999% Based On Trace Metals Analysis	Meets Requirements
Trace Metal Analysis	≤ 15.0 ppm

Specification: PRD.1.ZQ5.10000066095

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.