

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Lithium bromide solution - 54 wt. % in H₂O

Product Number:

411515

MDL:

MFCD00011077

Formula:

BrLi

Formula Weight:

86.85 g/mol

LiBr

TEST	Specification
------	---------------

Appearance (Clarity)

Clear

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Titration by AgNO₃

49.0 - 59.0 %

wt.% Solution As Lithium Bromide

ICP Major Analysis

Confirmed

Confirms Lithium Component

Specification: PRD.0.ZQ5.10000089815

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.