

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Haloperidol - powder

Product Number: H1512
CAS Number: 52-86-8
MDL: MFCD00051423
Formula: C₂₁H₂₃ClFNO₂
Formula Weight: 375.86 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Carbon	66.1 - 68.9 %
Nitrogen	3.4 - 4.0 %
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear
30 mg/ml, CHCl ₃	
Purity (TLC)	≥ 98 %
Recommended Retest Period	-----
3 years	

Specification: PRD.0.ZQ5.10000042870

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.