

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

(R)-(-)-6,6'-Dibromo-1,1'-bi-2-naphthol - 98%

Product Number:

482617

CAS Number:

65283-60-5

MDL:


MFCD00798290

Formula:

C₂₀H₁₂Br₂O₂

Formula Weight:

444.12 g/mol

**TEST****Specification**

Appearance (Color)

White to Yellow

Appearance (Form)

Powder or Crystals

Infrared spectrum

Conforms to Structure

Optical Rotation

-52.0 - -46.0 deg

c = 1.8%, THF

Purity (HPLC)

≥ 97.5 %

Specification: PRD.0.ZQ5.10000044858

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.