

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Glass beads, unwashed – 150–212 µm (70–100 U.S. sieve)

Product Number: **G9018**

TEST	Specification
Appearance (Color)	White to Light Blue
Appearance (Form)	Beads
Mesh Size	≥ 90 %
Within 150 - 212 microns	
Specification: PRD.0.ZQ5.10000045787	

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.