


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1,3,6,8-Tetrabromopyrene - 97%

Product Number: 717274
CAS Number: 128-63-2
Formula: C₁₆H₆Br₄
Formula Weight: 517.83 g/mol



TEST

Specification

Appearance (Color)	Off White to Yellow and Faint Beige to Beige and Faint Green to Green
Appearance (Form)	Conforms to Requirements
Powder or Solid	
Infrared Spectrum	Conforms to Structure
Titration by AgNO ₃	96.5 - 103.5 %
After Oxygen Combustion	

Specification: PRD.1.ZQ5.10000037081

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.