


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Nickel(II) acetylacetonate - 95%

Product Number: **283657**
CAS Number: 3264-82-2
MDL: MFCD00000024
Formula: C₁₀H₁₄NiO₄
Formula Weight: 256.91 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
Green	
Appearance (Form)	Powder or Crystals
Infrared Spectrum	Conforms to Structure
Titration with EDTA	21.5 - 24.2 %
% Ni	

Specification: PRD.3.ZQ5.10000026494

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.