

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Lycopene - $\geq 98\%$ (HPLC), from tomato

Product Number: **SMB00706**



TEST	Specification
Appearance (Form)	Powder
Appearance (Colour)	Red to Very Dark Red and Red-Brown and Very Dark Red-Brown
Solubility (Solvent)	Chloroform
Solubility (Conc)	0.95 - 1.05 mg/ml
Solubility (Turbidity)	Clear
Solubility (Color)	Orange to Red
Proton NMR spectrum	Conforms to Structure
Purity (HPLC)	$\geq 98.0\%$
Storage Conditions	Protect From Light
Recommended Retest Period	-----
2 Years	

Specification: PRE.0.ZQ5.10000018281

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.