

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Thiazole-4-carboxaldehyde - 90%

Product Number: **681105**
CAS Number: 3364-80-5

Formula: C₄H₃NOS
Formula Weight: 113.14 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)

White to Light Yellow and Faint Orange
to Light Orange and Faint Beige to Light Beige and Faint Brown
to Light Brown

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Titration by Hydroxylamine

88.0 - 112.0 %

Specification: PRD.0.ZQ5.10000045944

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.